


DNA/RNA extactions products

Xtralyz®

The DNA extraction chip. Our consumable that engages a large volume of sample before extracting DNA and RNA from all organisms trapped in it, even with complex matrices.

Provides efficient mechanical lysis without sample dilution or hazardous chemicals.





The Grinder

To perform the lysis and manage the fluids, the Grinder is a portable system ensuring an easy use of the Xtralyz[®].

Designed to ensure the proper lysis execution, it is also equipped with a waste container.

Direct'HD

An automatic system to process a large number of Xtralyz[®]. Automatically performs the lysis and wash step. Dedicated to users with a high daily test flow.

